

TalkBox Simulator User Manual

This system is very straight forward and easy to use. It has many features that allow the user to do many things.

When first looking at the setup, one can see that there are **three** buttons.

1. The **Recording Panel** button: This is where audio can be recorded. The purpose of this button is to record audio that can be used when pressing one of the buttons. After recording the user must put the audio file in a directory.
2. The **Setup Butons** button: This is where one can setup the orientation of the six buttons. In this setting/configuration panel the user can setup which audio goes to which button. They can also search for an existing audio by typing the word or they can use emoji's, i.e. smile (☺) or toilet (🚽). The user can also name the button and change its color. This is also where they can play he audio.
3. The **Create New Configuration Directory** button: This is where the user can create a new directoy for the configueration/setup. For example, if one likes their current configuraion execept for some things, they can make changes and close it. The new configuration will be saved if a new directory.

A "back" button was also created so that user can go back to the previous panel.

On the next page are a few screen shots of the system and its buttons:

TalkBox Configurer

Recording

Back to main menu

Filename of recording:rec.wav

Record

Stop Recording

TalkBox Configurer

Recording

Back to main menu

Filename of recording:rec.wav

Record

Stop Recording



